















# TARGON® INTERLOCKING NAIL SYSTEM »FEMUR«



| FEMUR<br>Universal      |            |                       |
|-------------------------|------------|-----------------------|
|                         | ø 10-11 mm | slotted<br>ø 12-15 mm |
| adapter                 | <b>A</b>   | <b>B</b>              |
| adapter bolt            | <b>A</b>   | <b>B</b>              |
| locking screw, proximal | ø 6,0 mm   | ø 6,0 mm              |
| drill, proximal         | ø 4,5 mm   | ø 4,5 mm              |
| locking screw, distal   | ø 5,0 mm   | ø 6,0 mm              |
| drill, distal           | ø 3,5 mm   | ø 4,5 mm              |
| sealing screw           | ø 8,0 mm   | ø 10,0 mm             |
| extractor bolt          | ø 8,0 mm   | ø 10,0 mm             |



|                            | TIBIA Universal |            | TIBIA Titanium |
|----------------------------|-----------------|------------|----------------|
|                            | slotted         |            | sdid           |
|                            | ø 10-11 mm      | ø 12-14 mm | ø 8-9 mm       |
| adapter                    | A               | B          | A              |
| adapter bolt               | A               | B          | A              |
| locking screw, dist./prox. | ø 5,0 mm        | ø 5,0 mm   | ø 4,5 mm       |
| drill                      | ø 3,5 mm        | ø 3,5 mm   | ø 3,5 mm       |
| sealing screw              | ø 8,0 mm        | ø 10,0 mm  | ø 8,0 mm       |
| extractor bolt             | ø 8,0 mm        | ø 10,0 mm  | ø 8,0 mm       |